


REGENETEN SOLUTION



Changing the
course of rotator
cuff disease

 **smith&nephew**
REGENETEN[®]
Bioinductive Implant

Supporting healthcare professionals

What is a rotator cuff tear?

The rotator cuff is a group of muscles and tendons that hold the shoulder joint in place and allow you to move your arm and shoulder. Rotator cuff tears are most often caused by chronic wear and tear with degeneration of the tendon and the prevalence of tears increases with age.¹ Among the most common musculoskeletal conditions, rotator cuff tendon tears are therapeutically challenging for orthopedic surgeons.

Risk Factors

Increasing age may represent a poor healing environment with reduced biomechanical properties of the tendon.²


Key Symptoms^{3,4}

- Pain at rest and at night, particularly if lying on the affected shoulder
- Pain when lifting and lowering your arm or with specific movements
- Weakness when lifting or rotating your arm
- Crepitus or crackling sensation when moving your shoulder in certain positions


Disease progression

- Rotator cuff disease is a progressive disease that worsens if left untreated:


Severe Tendinosis



Partial-Thickness Tear




Full-Thickness Tear



- Untreated rotator cuff tendinosis can progress to a partial – or full – thickness tear⁵
- Small tears tend to grow in size and severity, eventually requiring surgery⁶⁻⁸
- Up to 80% of partial-thickness tears increase in size within 2 years⁹

>50%


of asymptomatic rotator cuff tears progress in size within three years.¹³




The REGENETEN[◇] Solution

Traditional approaches focus only on biomechanical repair, but do not address the underlying biology of the tendon. The Smith & Nephew Rotator Cuff System is a new technology for improving the treatment of rotator cuff disease. It includes a bio-inductive implant that addresses both the biomechanics and biology required to heal a rotator cuff tear.¹²

Recovery: Current Standard of Care vs REGENETEN Solution



17 million people with rotator cuff disease in the U.S. at risk of disability.¹³





Key Outcomes

- Early intervention to reverse tear progression^{6, 7, 10}
- Preserves healthy tissue^{6, 10}
- Consistent healing of the tear^{6, 10}
- Fast return to normal activity¹³
- Potentially reduces re-tears¹⁰
- Helps to restore the rotator cuff footprint¹⁰

Bioinductive implant

This proprietary implant from Smith & Nephew stimulates the body's natural healing response to support new tendon growth and disrupt disease progression.⁶



The REGENETEN Bioinductive Implant gradually absorbs within 6 months and leaves a layer of new, tendon-like tissue to biologically augment the existing tendon.¹²

Sports Medicine
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810

www.smith-nephew.com
T +978 749 1000
US Customer Service:
+1 800 343 5717

*Trademark of Smith & Nephew.
©2019 Smith & Nephew.
All rights reserved. Printed in USA.
15920 V1 02/19

Supporting healthcare professionals for over 150 years

References

1. OrthoInfo (AAOS). Diseases & Conditions: Rotator cuff tears. 2017. <https://orthoinfo.aaos.org/en/diseases--conditions/rotator-cuff-tears/>. Accessed August 2018.
2. Mall N, Tanaka M, Choi L, Paletta Jr G. Factors Affecting Rotator Cuff Healing. *J Bone & Joint Surgery*. 2014; 96:778-88.
3. Mayo Clinic. Disease & conditions: Rotator cuff injury. 2017. <https://www.mayoclinic.org/diseases-conditions/rotator-cuff-injury/symptoms-causes/syc-20350225>. Accessed September 2018.
4. OrthoInfo (AAOS). Treatment: Rotator cuff tears: Surgical treatment options. 2017. <https://orthoinfo.aaos.org/en/treatment/rotator-cuff-tears-surgical-treatment-options/>. Accessed September 2018.
5. Keener JD, Galatz LM, Teefey SA et al. A prospective evaluation of survivorship of asymptomatic degenerative rotator cuff tears. *J Bone Joint Surg Am*. 2015;97:89-98.
6. Schlegel TF, Abrams JS, Bushnell BD, Brock JL, Ho CP. Radiologic and clinical evaluation of a bioabsorbable collagen implant to treat partial-thickness tears: a prospective multicenter study. *J Shoulder Elbow Surg*. 2018;27(2):242-251. doi: <http://dx.doi.org/10.1016/j.jse.2017.08.023>.
7. Bokor DJ, Sonnabend D, Deady L et al. Evidence of healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a collagen implant: a 2-year MRI follow-up. *MLTJ*. 2016;6(1):16-25.
8. Washburn R, Anderson TM, Tokish JM. Arthroscopic rotator cuff augmentation: Surgical technique using bovine collagen bioinductive implant. *Arthroscopy Techniques*. 2017;6(2):e297-e301.
9. Yamanaka K and Matsumoto T. The joint side tear of the rotator cuff: A followup study by arthrography. *Clin Orthop Relat Res*. 1994;(304):68-73.
10. Bokor DJ, Sonnabend D, Deady L et al. Preliminary investigation of a biological augmentation of rotator cuff repairs using a collagen implant: a 2-year MRI follow-up. *MLTJ*. 2015;5(3):144-150.
11. REBUILD Registry, data on file N=367.
12. Van Kampen C, Arnoczky, S, Parks, P, et al. Tissue-engineered augmentation of a rotator cuff tendon using a reconstituted collagen scaffold: A histological evaluation in sheep. *MLTJ*. 2013;3(3):229-235.
13. Yamaguchi K, Konstantino D, Middleton W et al. The Demographic and Morphological Features of Rotator Cuff Disease: A Comparison of Asymptomatic and Symptomatic Shoulders. *J Bone & Joint Surgery*. 2006 Aug; 88(8): 1699-1704.